

For life science research only.

Not for use in diagnostic procedures.

FOR *IN VITRO* USE ONLY.

Diffinity
Genomics®



FASTER, SIMPLER NUCLEIC ACID
ISOLATION & ANALYSIS

KANEKA

Diffinity DNA RapidExtract™ Kit

For rapid cell lysis for PCR-ready DNA

Version 2.0

Easy to use reagent kit for extracting DNA from sample material prior to PCR.

Part No. REXD-50/250	Kit for 50/250 reactions
----------------------	--------------------------

Store the kit at room temperature (+15 to +25 °C) protected from light.



-  This product contains alkaline reagents. Please observe general laboratory precautions and wear personal protective equipment.
-  Follow federal, state, and local regulations, and ensure safety when using or disposing of this product.

1. Diffinity DNA RapidExtract Kit

Number of Tests

The kit is designed for:

- 250 extraction reactions of blood samples
- 50 extraction reactions of tissue, cultured cells, bacteria, stool, or plant (leaf and seed) samples.

Kit Contents

Vial	Label	Contents
A	Solution A	5 ml Lysis Reagent
B	Solution B	0.7 ml Stop Solution

Storage and Stability

Store the Kit at +15 to +25 °C protected from light until the expiration date printed on the label.

Application

This kit is designed to easily obtain PCR-ready DNA extract from blood, animal tissues, cultured cells, plants, stool, and other sample types.

Assay Time

Procedure	Time
Combine sample with Solution A and mix.	1 minute
Incubate at 98 °C	8 minutes
Cool to Room temperature	2 minutes
Add Solution B and mix	1 minute
Total Time for Extraction	12 minutes

2. How this Product is Used

Leaf or Seed

1. Put a 5-8 mm² piece of plant material to a PCR tube. Add 100 µl of Solution A to the tube and mix well by pipetting.

* The sample size may vary depending on species and condition of the sample.

2. Homogenize the sample in the mixture, such as by crushing it with a pipette tip.

3. Incubate the PCR tube at 98°C for 8 minutes in a heat block.

* Secure the cap with a cap lock prior to incubation.

4. Cool the tube to room temperature.

* Make sure the tube has cooled down before opening the cap.

5. Add 14 µl of Solution B and mix well.

* Mix the extract from Step 5 before each use.

* If excessive amounts of solid residue are found in the extract, centrifuge the tube at 5000 rpm for 5 minutes at 4°C.

Use 1 to 5 µl of extracted template DNA in 50 µl PCR reaction.

Blood

1. Put 1 to 3 µl of blood to a PCR tube. Add 20 µl of the Solution A and mix well by pipetting.

2. Incubate the PCR tube at 98°C for 8 minutes in a heat block.

* Secure the cap with a cap lock prior to incubation.

4. Cool the tube to room temperature.

* Make sure the tube has cooled down before opening the cap.

5. Add 3 µl of Solution B and mix well.

* Mix the extract from Step 5 before each use.

* If excessive amounts of solid residue are found in the extract, centrifuge the tube at 5000 rpm for 5 minutes at 4°C.

Use 1 to 5 µl of extracted template DNA in 50 µl PCR reaction.

Tissue (mouse tail)

1. Put a 5-8 mm fragment of mouse tail into a PCR tube. Add 100 µl of Solution A and mix well by pipetting.

* The amount of starting material may vary depending on species and condition of the sample.

2. Homogenize the sample in the mixture, such as by crushing it with a pipette tip.

3. Incubate the PCR tube at 98°C for 8 minutes in a heat block.

* Secure the cap with a cap lock prior to incubation.

4. Cool the tube to room temperature.

* Make sure the tube has cooled down before opening the cap.

5. Add 14 µl of Solution B and mix well.

* Mix the extract from Step 5 before each use.

* If excessive amounts of solid residue are found in the extract, centrifuge the tube at 5000 rpm for 5 minutes at 4°C.

Use 1 to 5 µl of extracted template DNA in 50 µl PCR reaction.

Cultured cells or Microorganisms

1. Prepare a suspension of $10^3 - 10^5$ cells in a PCR tube.
 2. Centrifuge the tube and remove the supernatant. Add 100 μ l of Solution A to the cell pellet and mix well by pipetting.
 3. Incubate the PCR tube at 98°C for 8 minutes in a heat block.
 - * Secure the cap with a cap lock prior to incubation.
 4. Cool the tube to room temperature.
 - * Make sure the tube has cooled down before opening the cap.
 5. Add 14 μ l of Solution B and mix well.
 - * Mix the extract from Step 5 before each use.
 - * If excessive amounts of solid residue are found in the extract, centrifuge the tube at 5000 rpm for 5 minutes at 4°C.
- Use 1 to 5 μ l of extracted template DNA in 50 μ l PCR reaction.

Stool

1. Prepare suspension of stool sample in 100 μ l sterilized water. Put 10 μ l of the suspension into a PCR tube. Add 100 μ l of Solution A to the suspension and mix well by pipetting.
 - * The starting amount of stool may vary depending on species and condition of the sample.
 3. Incubate the PCR tube at 98°C for 8 minutes in a heat block.
 - * Secure the cap with a cap lock prior to incubation.
 4. Cool the tube to room temperature.
 - * Make sure the tube has cooled down before opening the cap.
 5. Add 14 μ l of Solution B and mix well.
 - * Mix the extract from Step 5 before each use.
 - * If excessive amounts of solid residue are found in the extract, centrifuge the tube at 5000 rpm for 5 minutes at 4°C.
- Use 1 to 5 μ l of extracted template DNA in 50 μ l PCR reaction.

Diffinity DNA RapidExtract Kit - Manual Version 2 (December 2015)

Diffinity Genomics | Chiral Technologies, Inc.

800 North Five Points Road

West Chester, PA 19380

+1 800-624-4725

www.diffinitygenomics.com

©2015 Diffinity Genomics | Chiral Technologies. All rights reserved. Diffinity, Diffinity Genomics, and RapidExtract are trademarks of Chiral Technologies, Inc.